

## RBC Lysis Buffer (10X)

**Catalog # / Size:** 420301 / 100 ml

**Preparation:** The pH of the 1X solution should fall within the range of pH 7.1-7.4. Adjust the pH if necessary. Warm the 1X solution to room temperature prior to use.

**Formulation:** This red blood cell (RBC) lysis buffer is supplied as a 10X solution and should be diluted to 1X in deionized water. 100 ml of 10X concentrate will yield a quantity of 1X solution that is sufficient to lyse 500 samples.

**Storage:** Store RBC lysis buffer between 2°C and 8°C.

## Applications:

**Applications:** FC, ICFC - *Quality tested*

**Application Notes:** **A. Lysis of Mouse Spleen RBCs:**

1. Harvest mouse spleen and prepare a single cell suspension.
2. Pellet the cells by centrifugation (350 x g); aspirate the supernatant.
3. Dilute the 10X Red Blood Cell Lysis Buffer to 1X working concentration with deionized water and resuspend the pellet in 5 ml of 1X Lysis Buffer.
4. Incubate on ice for 4-5 minutes with occasional shaking.
5. Stop the reaction by diluting the Lysis Buffer with 20-30 ml of 1X PBS.
6. Spin the cells (350 x g) and discard the supernatant.
7. Resuspend the pellet in the appropriate buffer (e.g., BioLegend Cell Staining Buffer Cat. No. 420201), wash 1X.
8. Count cells, adjust density, and proceed with cell staining procedures.

**B. Lysis of Human Peripheral Blood RBCs:**

1. Dilute the 10X RBC Lysis Buffer to 1X working concentration with deionized water. Warm the 1X solution to room temperature prior to use.
2. Add 2.0 ml of 1X RBC Lysis Buffer to each tube containing up to 100 µl of whole blood.
3. Gently vortex each tube immediately after adding the lysing solution. Incubate at room temperature, protected from light, for 10-15 minutes.
4. Centrifuge 350 x g for 5 minutes. Aspirate supernatant without disturbing pellet.
5. Resuspend the pellet in the appropriate buffer (e.g., BioLegend Cell Staining Buffer, Cat. No. 420201), wash 1X.
6. Resuspend and proceed with further procedures.

### Application References:

1. Adams RA, *et al.* 2007. *J. Exp. Med.* 204:571.
2. Feild-Corbett C, *et al.* 2009. *Stem Cells.* PubMed
3. de Resende MM, *et al.* 2010. *Physiol Genomics.* 42:437. PubMed
4. Roybal JD, *et al.* 2011. *Mol Cancer Res.* 9:25. PubMed
5. Boldin MP, *et al.* 2011. *J Exp Med.* PubMed
6. Kahlenberg JM, *et al.* 2011. *J. Immunol.* 187:6143. PubMed
7. Deiuliis JA, *et al.* 2012. *Am J Physiol Lung Cell Mol Physiol.* 302:399. PubMed
8. Wu CY, *et al.* 2012. *PLoS One.* 7:e42365. PubMed
9. Hsia BJ, *et al.* 2012. *J Allergy Clin Immunol.* 130:205. PubMed
10. Rico C, *et al.* 2012. *Carcinogenesis.* 33:2283. PubMed
11. Jin J, *et al.* 2012. *Arterioscler Thromb Vasc Biol.* 32:2901. PubMed
12. Saraiva AM, 2013. *Hum Immunol.* 74:207. PubMed
13. Koon HW, *et al.* 2013. *J Vis exp.* 68:4208. PubMed
14. Ma Dy, *et al.* 2013. *J Virol.* 87:3361. PubMed
15. Han W, 2013. *J. Immunol.* 190:4786. PubMed
16. Garcia-Acros I, *et al.* 2013. *J Biol Chem.* 288:14046. PubMed
17. Zhao JL, *et al.* 2013. *Elife.* 21:537. PubMed
18. Guo W, *et al.* 2013. *J Gerontol A Biol Sci Med Sci.* PubMed
19. Yamasaki S, *et al.* 2013. *Clin Vaccine immunol.* 20:1508. PubMed
20. Kraft P, *et al.* 2013. *Stroke.* 44:3202. PubMed
21. Liu X, *et al.* 2013. *Invest Ophthalmol Vis Sci.* 54:7386. PubMed
22. Fullerton AM, *et al.* 2013. *Toxicol Sci.* 136:72. PubMed
23. Cekic C, *et al.* 2013. *J Exp Med.* 210:2693. PubMed
24. Varga T, *et al.* 2013. *J. Immunol.* 191:5695. PubMed
25. Ma Y, *et al.* 2014. *Cancer Res.* 74:436. PubMed
26. Chang YJ, *et al.* 2014. *Cytotherapy.* 16:560. PubMed
27. Kim MJ, *et al.* 2014. *Mol Pharmacol.* 86:284. PubMed
28. Cekic C, *et al.* 2014. *Cancer Res.* 74:7250. PubMed

**Description:** Red Blood Cell (RBC) Lysis Buffer has been designed, formulated, and tested to ensure optimal lysis of RBCs in single cell suspensions with minimal effects on leukocytes. RBC Lysis Buffer is supplied as a 10X solution containing ammonium chloride, potassium carbonate, and EDTA, and should be diluted in deionized water prior to use. Nucleated RBCs are not effectively lysed with ammonium chloride.



For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.



\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biollegend.com/ordering#license](http://www.biollegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

**Related Products:Product**  
Cell Staining Buffer

**Clone**

**Application**  
FC, ICC, ICFC



For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.



\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 9727 Pacific Heights Blvd, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587