



MojoSort[™] Human anti-APC Nanobeads



A single cell suspension from peripheral blood mononuclear cells (PBMCs) was prepared and stained with TCR V δ 2 (clone B6) APC. APC-positive cells were then isolated using MojoSort[™] Human anti-APC Nanobeads. Dead cells were excluded by 7AAD.

MojoSort[™] Human anti-APC Nanobeads



A single cell suspension from peripheral blood mononuclear cells (PBMCs) was prepared and stained with CD3 (clone UCHT1) APC. APC-positive cells were then isolated using MojoSort[™] Human anti-APC Nanobeads. Dead cells were excluded by 7AAD.