

## Brefeldin A Solution (1,000X)

<b>Catalog# / Size</b>	420601 / 1 mL
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	BFA
<b>Description</b>	Brefeldin A (BFA) is a protein transport inhibitor commonly used to enhance intracellular cytokine staining signals by blocking transport processes during cell activation. Especially useful for the intracellular staining of cytokines, BFA leads to the accumulation of most cytokines at the Golgi Complex/Endoplasmic Reticulum (see Jung, <i>et al.</i> , 1993). Optimal conditions for use are cell type and time-dependent. Typically, protein transport inhibitors are included during <i>in vitro</i> cell activation cultures for 4-24 hours prior to harvest (see references below for additional information). Brefeldin A Solution is supplied as a 1,000X solution, which should be diluted to 1X in cell culture medium.

### Product Details

---

<b>Formulation</b>	Brefeldin A Solution is supplied as a 1000X in DMSO, which should be diluted to 1X in cell culture medium.
<b>Concentration</b>	5.0 mg/ml
<b>Storage &amp; Handling</b>	Store the Brefeldin A Solution (1,000X) between 2°C and 8°C. Note: DMSO freezes at this temperature.
<b>Application</b>	<a href="#">ICFC - Quality tested</a>
<b>Recommended Usage</b>	Dilute the 1000X solution to 1X in the tissue culture medium. It is recommended that cells are cultured with brefeldin A for ≤ 24 hours, as this can become toxic for cell viability.
<b>Application Notes</b>	Completely thaw the solution before use, DMSO freezes at 4°C. We recommend aliquotting the reagent into smaller volumes to avoid repeated freeze-thawing.

### Application References

- O'Sullivan BJ, *et al.* 2006. *J. Immunol.* 176:7278.
- Fuse S, *et al.* 2007. *J. Immunol.* 178:5227.
- Kang YJ, *et al.* 2007. *Nature Immunol.* 8:601.
- Redfern CH, *et al.* 2006. *J. Clin Oncol.* 19:3107. [PubMed](#)
- Smithy MJ, *et al.* 2008. *J. Immunol.* 180:3406. [PubMed](#)
- Elzey BD, *et al.* 2008. *Blood* 111:38684. [PubMed](#)
- Nguyen KD, *et al.* 2008. *J. Immunol.* 181:5386. [PubMed](#)
- Miner MD, *et al.* 2008. *Microb Pathog.* 45:273. [PubMed](#)
- Goodridge HS, *et al.* 2009. *J. Immunol.* 182:1146. [PubMed](#)
- Markey KA, *et al.* 2009. *Blood* 113:5644. [PubMed](#)
- Wilcox RA, *et al.* 2009. *Blood* [PubMed](#)
- Zenaro E, *et al.* 2009. *J. Leukoc Bio.* [PubMed](#)
- Rozkova D, *et al.* 2009. *Clin Immunol.* 131:1. [PubMed](#)
- Butcher JP, *et al.* 2010. *Clin Cancer Res.* 16:2065.
- Boisvert M, *et al.* 2010. *J. Virol.* 84:7782. [PubMed](#)
- Liu WC, *et al.* 2010. *J. Virol.* 84:12011.15. Tawara I, *et al.* 2011. *Clin Cancer Res.* 17:77. [PubMed](#)
- Iwata Y, *et al.* 2011. *Blood.* 117:530. [PubMed](#)
- Schmidt NW, *et al.* 2011. *J. Immunol.* 186:3836. [PubMed](#)
- Langkjaer A, *et al.* 2012. *Clin Immunol.* 145:1. [PubMed](#)
- Molinero LL, *et al.* 2012. *PNAS.* 109:18529. [PubMed](#)
- Li HW, *et al.* 2012. *J. Immunol.* 189:5572. [PubMed](#)
- Wang JM, *et al.* 2013. *J. Virol.* 87:4372. [PubMed](#)
- Fukuoka A, *et al.* 2013. *Int Immunol.* 25:373. [PubMed](#)
- Wang Jm, *et al.* 2013. *J Virol.* 87:11626. [PubMed](#)
- Lee TH, *et al.* 2013. *J Biol Chem.* 288:32873. [PubMed](#)
- Shi L, *et al.* 2014. *J. Immunol.* 192:649. [PubMed](#)
- Bulati M, *et al.* 2014. *Exp Gerontol.* [PubMed](#)
- Blankenhaus B, *et al.* 2014. *PLoS Pathog.* 10:1003913. [PubMed](#)
- Tuna H, *et al.* 2014. *J Leukoc Biol.* 95:471. [PubMed](#)
- Miller ML, *et al.* 2014. *PNAS.* 111:7397. [PubMed](#)
- Wang Y, *et al.* 2014. *Sci Rep.* 4:5346. [PubMed](#)

32. Zhou GL, *et al.* 2014. *J Cell Sci.* 127:3916. [PubMed](#)
33. Rist MJ, *et al.* 2015. *J Virol.* 89:703. [PubMed](#)

## Product Citations

1. Rossato C, *et al.* 2017. *Biologicals.* 10.1016/j.biologicals.2017.06.007. [PubMed](#)
2. Ruland C, *et al.* 2017. *Brain Behav Immun.* 10.1016/j.bbi.2017.06.010. [PubMed](#)
3. Zischke J, *et al.* 2017. *PLoS Pathogens.* 13(6):e1006454. [PubMed](#)
4. Ries M, *et al.* 2017. *PLoS Pathog.* 10.1371/journal.ppat.1006506. [PubMed](#)
5. Martins M, *et al.* 2017. *PLoS Pathog.* 10.1371/journal.ppat.1006529. [PubMed](#)
6. Dayoub O, *et al.* 2017. *Sci Rep.* 10.1038/s41598-017-09933-4. [PubMed](#)
7. Orlov M, *et al.* 2017. *PLoS One.* 10.1371/journal.pone.018296. [PubMed](#)
8. Li L, *et al.* 2017. *Sci Rep.* 10.1038/s41598-017-14000-z. [PubMed](#)
9. Ahmadi H, *et al.* 2017. *Inflammopharmacology.* 10.1007/s10787-017-0410-8. [PubMed](#)
10. Maschmeyer P, *et al.* 2017. *J Autoimmun.* 10.1016/j.jaut.2017.11.005. [PubMed](#)
11. Landi A, *et al.* 2017. *Vaccine.* 10.1016/j.vaccine.2017.10.072. [PubMed](#)
12. Zhu S, *et al.* 2017. *Exp Cell Res.* 10.1016/j.yexcr.2017.10.013. [PubMed](#)
13. Orozco SL, *et al.* 2019. *Cell Rep.* 28:2275. [PubMed](#)
14. Yuan X, *et al.* 2017. *Elife.* 6:e29540. [PubMed](#)
15. Lu G, *et al.* 2017. *PLoS One.* 12:e0178352. [PubMed](#)
16. Gebremeskel S, *et al.* 2017. *Cancer Immunol Res.* 0.9625. [PubMed](#)
17. Li J, *et al.* 2018. *Redox Biol.* 0.976388889. [PubMed](#)
18. Cao Y, *et al.* 2018. *Nat Commun.* 2.940972222. [PubMed](#)
19. Emmer BT, *et al.* 2018. *Elife.* 7:e38839. [PubMed](#)
20. Guo C, *et al.* 2018. *Immunity.* 49:842. [PubMed](#)
21. Wang R, *et al.* 2018. *Oncoimmunology.* 7:e1461303. [PubMed](#)
22. Carissimo G, *et al.* 2019. *Life Sci Alliance.* 2:e201900298. [PubMed](#)
23. Salazar V, *et al.* 2019. *Cell Rep.* 26:1585. [PubMed](#)
24. Ngwenyama N, *et al.* 2019. *JCI Insight.* 4:e125527. [PubMed](#)
25. Yamamoto A, *et al.* 2019. *J Clin Invest.* 129:3640. [PubMed](#)
26. Piper CJM, *et al.* 2020. *Cell Reports.* 29(7):1878-1892.e7. [PubMed](#)
27. Raju S, *et al.* 2020. *Cell Reports.* 29(9):2556-2564.e3. [PubMed](#)
28. Rana A, *et al.* 2020. *Front Immunol.* 11:536. [PubMed](#)
29. Good Z, *et al.* 2019. *Nat Biotechnol.* 37:259. [PubMed](#)
30. Gogokhia L, *et al.* 2019. *Cell Host Microbe.* 25:285. [PubMed](#)
31. Rosser EC, *et al.* 2020. *Cell Metabolism.* 31(4):837-851. [PubMed](#)
32. Hassan AO, *et al.* 2020. *Cell.* 183:169. [PubMed](#)
33. Sekine T, *et al.* 2020. *Cell.* 183:158. [PubMed](#)
34. Wang Y, *et al.* 2020. *Sci Adv.* 6:eaba5412. [PubMed](#)
35. Kalafati L, *et al.* 2020. *Cell.* 183:771. [PubMed](#)
36. Redfern C, *et al.* 2006. *J Clin Oncol.* 24:3107. [PubMed](#)
37. Nguyen K, *et al.* 2008. *J Immunol.* 181:5386. [PubMed](#)
38. Elzey B, *et al.* 2008. *Blood.* 111:3684. [PubMed](#)
39. Smithy M, *et al.* 2008. *J Immunol.* 180:3406. [PubMed](#)
40. Zenaro E, *et al.* 2009. *J Leukoc Biol.* 86:1393. [PubMed](#)
41. Wilcox R, *et al.* 2009. *Blood.* 114:2936. [PubMed](#)
42. Markey K, *et al.* 2009. *Blood.* 113:5644. [PubMed](#)
43. Goodridge H, *et al.* 2009. *J Immunol.* 182:1146. [PubMed](#)
44. Liu W, *et al.* 2010. *J Virol.* 84:12011. [PubMed](#)
45. Gnerlich J, *et al.* 2010. *J Immunol.* 185:4063. [PubMed](#)
46. Boisvert M, *et al.* 2010. *J Virol.* 84:7782. [PubMed](#)
47. Butchar J, *et al.* 2010. *Clin Cancer Res.* 2.100694444. [PubMed](#)
48. Yu H, *et al.* 2011. *J Leukoc Biol.* 90:205. [PubMed](#)
49. Harty N 2011. *J Immunol.* 186:3836. [PubMed](#)
50. Iwata Y, *et al.* 2011. *Blood.* 117:530. [PubMed](#)
51. Tawara I, *et al.* 2011. *Clin Cancer Res.* 0.761805556. [PubMed](#)
52. Langkjær A, *et al.* 2012. *Clin Immunol.* 145:1. [PubMed](#)
53. Molinero L, *et al.* 2012. *Proc Natl Acad Sci U S A.* 109:18529. [PubMed](#)
54. Li H, *et al.* 2012. *J Immunol.* 189:5572. [PubMed](#)
55. Ji L, *et al.* 2012. *PLoS One.* 7:e50909. [PubMed](#)
56. Fukuoka A, *et al.* 2013. *Int Immunol.* 25:373. [PubMed](#)
57. Wang J, *et al.* 2013. *J Virol.* 87:4372. [PubMed](#)
58. Rozková D, *et al.* 2009. *Clin Immunol.* 131:1. [PubMed](#)
59. Miner M, *et al.* 2008. *Microb Pathog.* 45:273. [PubMed](#)
60. Wang J, *et al.* 2013. *J Virol.* 87:11262. [PubMed](#)
61. Lee T, *et al.* 2013. *J Biol Chem.* 288:32873. [PubMed](#)
62. Shi L, *et al.* 2014. *J Immunol.* 192:649. [PubMed](#)
63. Bulati M, *et al.* 2014. *Exp Gerontol.* 54:123. [PubMed](#)
64. Blankenhaus B, *et al.* 2014. *PLoS Pathog.* 10:1003913. [PubMed](#)
65. Tuna H, *et al.* 2014. *J Leukoc Biol.* 95:471. [PubMed](#)
66. Miller M, *et al.* 2014. *Proc Natl Acad Sci U S A.* 111:7397. [PubMed](#)
67. Wang Y, *et al.* 2014. *Sci Rep.* 4:5346. [PubMed](#)
68. Zhou G, *et al.* 2014. *J Cell Sci.* 127:3916. [PubMed](#)
69. Rist M, *et al.* 2015. *J Virol.* 89:703. [PubMed](#)
70. Vigano S, *et al.* 2015. *J Virol.* 89: 7829-7840. [PubMed](#)
71. Zhang W, *et al.* 2015. *Infect Immun.* 83: 3213-3223. [PubMed](#)
72. Starnbach X 2015. *J Immunol.* 195: 1665-1675. [PubMed](#)
73. Hoechst B, *et al.* 2016. *J Immunol.* 195:1517-1523. [PubMed](#)
74. Pablo-Bernal R, *et al.* 2016. *J Gerontol A Biol Sci Med Sci.* 71: 610 - 618. [PubMed](#)
75. Matsushita T, *et al.* 2015. *Rheumatology.* 101093/rumatology/kv331. [PubMed](#)

76. Gartlan K, *et al.* 2015. *Blood*. 126: 1609-1620. [PubMed](#)
77. Morse M, *et al.* 2016. *Lupus*. 25: 81 - 87. [PubMed](#)
78. Pablo-Bernal R, *et al.* 2015. *J Infect Dis*. 10.1093/infdis/jiv520. [PubMed](#)
79. Rau M, *et al.* 2016. *J Immunol*. 196: 97 - 105. [PubMed](#)
80. Krishna B, *et al.* 2016. *Sci Rep*. 6: 24674. [PubMed](#)
81. Garg B, *et al.* 2016. *Mutagenesis*. 10.1093/mutage/gew023. [PubMed](#)
82. Gonzalez-Perez G, *et al.* 2016. *J Immunol*. 196: 3768 - 3779. [PubMed](#)
83. Doll K, *et al.* 2016. *J Immunol*. 196: 4253 - 4262. [PubMed](#)
84. Rialdi A, *et al.* 2016. *Science*. 352: aad7993. [PubMed](#)
85. Satoh M, *et al.* 2016. *Sci Rep*. 6:28473. [PubMed](#)
86. Castro I, *et al.* 2016. *J Virol*. 90: 5280 - 5291. [PubMed](#)
87. Vacaflores A, *et al.* 2016. *PLoS One*. 11: 0157175. [PubMed](#)
88. Ma H, Wan S, Xia C 2016. *J Leukoc Biol*. 99: 1121 - 1129. [PubMed](#)
89. Ruhland M, *et al.* 2016. *Nat Commun*. 7:11762. [PubMed](#)
90. Nice T, *et al.* 2016. *PLoS One*. 12: 1005684. [PubMed](#)
91. Liu W, *et al.* 2016. *J Virol*. 90: 6085 - 6096. [PubMed](#)
92. Kim J, *et al.* 2016. *Nat Immunol*. 17:1159-1166. [PubMed](#)
93. Al-Attar A, *et al.* 2016. *Data Brief*. 8: 536-44. [PubMed](#)
94. Li B, Schmidt N 2016. *PLoS One*. 11: 0162427. [PubMed](#)
95. Krämer B, *et al.* 2016. *PLoS One*. 11: 0162068. [PubMed](#)
96. Ikebuchi R, *et al.* 2016. *Sci Rep*. 6:35002. [PubMed](#)
97. Gross C, *et al.* 2016. *Neurol Neuroimmunol Neuroinflamm*. 3: e289. [PubMed](#)
98. Ulland T, *et al.* 2016. *Nat Commun*. 7:13180. [PubMed](#)
99. Mackroth M, *et al.* 2016. *PLoS Pathog*. 12:e1005909. [PubMed](#)
100. Ciarlo E, *et al.* 2016. *Sci Rep*. 6:37944. [PubMed](#)
101. Pfeifle R, *et al.* 2016. *Nat Immunol*. 18:104-113. [PubMed](#)
102. Abshire C, *et al.* 2016. *PLoS Pathog*. 12:e1006088. [PubMed](#)
103. Fucikova J, *et al.* 2016. *Blood*. 128(26):3113-3124. [PubMed](#)
104. Griffiths K, *et al.* 2016. *Nat Commun*. 7:13894. [PubMed](#)
105. Kadivar M, *et al.* 2016. *J Immunol*. 197(12):4584-4592. [PubMed](#)
106. Beyaz S, *et al.* 2016. *Nat Immunol*. 18:184-195. [PubMed](#)
107. Just S, *et al.* 2016. *Sci Rep*. 6:39796. [PubMed](#)
108. Matisz C, *et al.* 2017. *Sci Rep*. 7:40631. [PubMed](#)
109. Douthwaite J, *et al.* 2017. *J Immunol*. 198(1):528-537. [PubMed](#)
110. Dáňová K, *et al.* 2017. *J Immunol*. 198(2):729-740. [PubMed](#)
111. Hradilova N, *et al.* 2017. *PLoS One*. 12(2):e0171539. [PubMed](#)
112. Schridde A, *et al.* 2017. *Mucosal Immunol*. 10.1038/mi.2016.142. [PubMed](#)
113. Wang W, *et al.* 2020. *STAR Protoc*. 1:100218. [PubMed](#)
114. Johnson AM, *et al.* 2020. *J Immunol*. 204:2295. [PubMed](#)
115. Shi Z, *et al.* 2021. *Ann Clin Transl Neurol*. 8:43. [PubMed](#)
116. Laczkó D, *et al.* 2020. *Immunity*. 53:724. [PubMed](#)
117. Buqué A, *et al.* 2020. *Nat Commun*. 3.110416667. [PubMed](#)
118. Festag J, *et al.* 2020. *Mol Ther Nucleic Acids*. 1.330555556. [PubMed](#)
119. Loeuillard E, *et al.* 2020. *J Clin Invest*. 130:5380. [PubMed](#)
120. Codo AC, *et al.* 2020. *Cell Metab*. 32:437. [PubMed](#)
121. Mohammadpour H, *et al.* 2020. *JCI Insight*. 5:00. [PubMed](#)
122. Yan H, *et al.* 2020. *ACS Appl Mater Interfaces*. 12:19324. [PubMed](#)
123. Tang B, *et al.* 2020. *Clin Cancer Res*. 26:2216. [PubMed](#)

## Antigen Details

---

### Antigen References

1. Current Protocols in Immunology (John Wiley & Sons New York) Unit 6.24 *Detection of Intracellular Cytokines by Flow Cytometry* (Barbara Foster and Calman Prussin NIAID NIH Bethesda MD).
2. Sander B, *et al.* 1991. *Immunol. Rev*. 119:65.
3. Sander B, *et al.* 1993. *J. Immunol. Meth*. 166:201.
4. Prussin C, *et al.* 1995. *J. Immunol. Meth*. 188:117.
5. Jung T, *et al.* 1993. *J. Immunol. Meth*. 159:197.

### Gene ID

NA

## Related Protocols

---

[Intracellular Cytokine Staining Protocol - Video](#)

[Intracellular Flow Cytometry Staining Protocol](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587