


## **Purified anti-Connexin 43, 360-382 Antibody**



IHC staining of purified anti-Connexin 43, 360-382 antibody (clone P2C4) on formalin-fixed, paraffin-embedded normal human brain tissue. Following antigen retrieval using Sodium Citrate H.I.E.R., the tissue was incubated with the primary antibody at 5 µg/ml overnight at 4°C. BioLegend's Ultra-Streptavidin (USA) HRP kit was used for detection followed by hematoxylin counterstaining, according to the protocol provided.