


Purified anti-Keratin 15 Antibody



IHC staining of anti-Keratin 15 antibody (Poly18339) on formalin-fixed paraffin-embedded human skin tissue. Following antigen retrieval using Sodium Citrate, the tissue was incubated with the primary antibody at 1 $\mu\text{g}/\text{ml}$ for 60 minutes at room temperature. After incubation with biotinylated anti-chicken antibody, the tissue was labeled with BioLegend's HRP labeling reagent followed by addition of AEC Chromogen Concentrate. The tissue was then counterstained with hematoxylin, according to the protocol provided.