


Purified anti-human IgD Antibody



SeqIF[™] (sequential immunofluorescence) staining on COMET[™] of Purified anti-IgD (clone W18340A, yellow) on formalin-fixed paraffin-embedded human tonsil tissue at 5 $\mu\text{g}/\text{mL}$. Alexa Fluor[™] Plus 555 Goat anti-Rat IgG antibody (Lunaphore, Cat. No. DR555RT) was used as a secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.