

LEGENDplex™ Mouse IL-6 Capture Bead B4, 13X

Catalog# / Size	740159 / 100 tests (270 µL)
Regulatory Status	RUO
Description	LEGENDplex™ Mouse Inflammation Mix and Match Panel Capture Beads are intended for use with the following reagents: 740371 (LEGENDplex™ Mouse Inflammation Panel Standard) 740165 (LEGENDplex™ Mouse Inflammation Panel Detection Antibodies) 740373 (LEGENDplex™ Buffer Set B)

Product Details

Verified Reactivity	Mouse
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Product Citations	<ol style="list-style-type: none">1. Meek SM, <i>et al.</i> 2018. Pathogens. 7:22. PubMed2. Aguayo–Mazzucato C, <i>et al.</i> 2019. Cell Metab. 30:129. PubMed
--------------------------	---

Antigen Details

Biology Area	Stem Cells, Immunology, Innate Immunity
Molecular Family	Cytokines/Chemokines
Gene ID	16193

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587

