

Purified anti-mouse I-A^k (A β ^k) Antibody

Catalog# / Size	109902 / 500 μ g
Clone	10-3.6
Regulatory Status	RUO
Other Names	MHC class II
Isotype	Mouse (CWB) IgG2a, κ
Description	The 10-3.6 antibody reacts with the β chain of the I-A ^k MHC class II alloantigen. This class II molecule is expressed on antigen presenting cells (including B cells) and a subset of T cells from H-2 ^k bearing mice and involved in antigen presentation to T cells expressing CD3/TCR and CD4 proteins. The 10-3.6 antibody cross-reacts with I-A ^{f,r,s} antigens and I-A ^{g7} of NOD mice; it does not react with other haplotypes (e.g., b, d, p, q).

Product Details

Verified Reactivity	Mouse
Antibody Type	Monoclonal
Host Species	Mouse
Immunogen	C3H mouse splenocytes
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Preparation	The antibody was purified by affinity chromatography.
Concentration	0.5 mg/ml
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C.
Application	FC - Quality tested IHC-F, IP - Reported in the literature, not verified in house
Recommended Usage	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis . For flow cytometric staining, the suggested use of this reagent is \leq 0.25 μ g per million cells in 100 μ l volume. It is recommended that the reagent be titrated for optimal performance for each application.
Application Notes	Additional reported applications (for the relevant formats) include immunoprecipitation ^{1,2} , protection against autoimmune IDDM ³ , <i>in vitro</i> blocking of antigen-specific MHC-restricted responses and <i>in vivo</i> inhibition of lymphoma growth ⁴ , and immunohistochemical staining ⁵ of acetone-fixed frozen sections. This clone does not cross react with other haplotypes (e.g. b, d, p, q).
Application References (PubMed link indicates BioLegend citation)	1. Landais D, <i>et al.</i> 1986. <i>J. Immunol.</i> 137:3002. (IP) 2. Kupinski JM, <i>et al.</i> 1983. <i>J. Immunol.</i> 130:2277. (IP) 3. Kappler JW, <i>et al.</i> 1981. <i>J. Exp. Med.</i> 153:1198. 4. Alisauskas RM, <i>et al.</i> 1986. <i>Immunopharmacology</i> 12:1. 5. Reis e Sousa and Germain 1999. <i>J. Immunol.</i> 162:6652. (IHC) 6. Yui MA, <i>et al.</i> 2010. <i>J. Immunol.</i> 185:284. PubMed 7. Gaudreau S, <i>et al.</i> 2007. <i>J. Immunol.</i> 179:3638. (FC)
Product Citations	1. Drake JR, <i>et al.</i> 2019. <i>Immunohorizons.</i> 3:28. PubMed
RRID	AB_313451 (BioLegend Cat. No. 109902)

Antigen Details

Structure	β chain of MHC class II
Distribution	B cell and activated T cells, APCs of H-2 ^b mice
Function	Antigen presentation
Ligand/Receptor	CD3/TCR, CD4
Cell Type	Antigen-presenting cells, B cells, T cells, Tregs
Biology Area	Immunology, Innate Immunity
Molecular Family	MHC Antigens
Antigen References	1. Watts C. 1997. <i>Ann. Rev. Immunol.</i> 15:821. 2. Pamer E, <i>et al.</i> 1998. <i>Ann. Rev. Immunol.</i> 16:323.
Gene ID	14960


Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

FITC anti-mouse I-A^k (A β ^k), PE anti-mouse I-A^k (A β ^k), Purified anti-mouse I-A^k (A β ^k)

Product Data



C3H mouse splenocytes stained with 10-3.6 purified, followed by anti-mouse IgG2a FITC

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587