

Alexa Fluor® 647 Rat IgG1, λ Isotype Ctrl Antibody

Catalog# / Size	401923 / 25 µg 401924 / 100 µg
Clone	G0114F7
Regulatory Status	RUO
Isotype	Rat IgG1, λ
Description	This antibody was chosen as an isotype control after screening on a variety of live mouse, rat, and human tissues.

Product Details

Antibody Type	Monoclonal
Host Species	Rat
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide
Preparation	The antibody was purified by affinity chromatography and conjugated with Alexa Fluor® 647 under optimal conditions.
Concentration	0.5 mg/mL
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	<p>Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest. Use our Concentration Lookup tool to find the exact concentrations of your lots of product.</p> <p>* Alexa Fluor® 647 has a maximum emission of 668 nm when it is excited at 633 nm / 635 nm.</p> <p>Alexa Fluor® and Pacific Blue™ are trademarks of Life Technologies Corporation.</p> <p>View full statement regarding label licenses</p>
Excitation Laser	Red Laser (633 nm)

Antigen Details

Gene ID	NA
----------------	----

Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

Purified Rat IgG1, λ Isotype Ctrl, APC Rat IgG1, λ Isotype Ctrl, PE Rat IgG1, λ Isotype Ctrl, PE/Cyanine7 Rat IgG1, λ Isotype Ctrl, PerCP/Cyanine5.5 Rat IgG1, λ Isotype Ctrl, Brilliant Violet 421™ Rat IgG1, λ Isotype Ctrl, FITC Rat IgG1, λ Isotype Ctrl, Ultra-LEAF™ Purified Rat IgG1, λ Isotype Ctrl, APC/Fire™ 750 Rat IgG1, λ Isotype Ctrl, TotalSeq™-A0237 Rat IgG1, λ Isotype Ctrl, TotalSeq™-

C0237 Rat IgG1, λ Isotype Ctrl, Alexa Fluor® 647 Rat IgG1, λ Isotype Ctrl

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587