

Brilliant Violet 605™ anti-mouse CD196 (CCR6) Antibody

Catalog# / Size	129819 / 50 µg
Clone	29-2L17
Regulatory Status	RUO
Other Names	CD196, BN-1, CKR-L3, DRY6, KY411, STRL22
Isotype	Armenian Hamster IgG
Description	CCR6 is a G-protein linked chemokine receptor which binds the chemokine CCL20/MIP-3a. It is expressed on B lymphocytes and certain subsets of dendritic cells and T cells. CCR6 is reported to be involved in mucosal immune response and lymphocyte migration and homeostasis.

Product Details

Verified Reactivity	Mouse
Antibody Type	Monoclonal
Host Species	Armenian Hamster
Immunogen	N-terminal peptide (aa. 2-38) of mouse CCR6 fused with GST
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA).
Preparation	The antibody was purified by affinity chromatography and conjugated with Brilliant Violet 605™ under optimal conditions.
Concentration	0.2 mg/ml
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	<p>Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤0.5 µg per million cells in 100 µl volume or 100 µl of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.</p> <p>Brilliant Violet 605™ excites at 405 nm and emits at 603 nm. The bandpass filter 610/20 nm is recommended for detection, although filter optimization may be required depending on other fluorophores used. Be sure to verify that your cytometer configuration and software setup are appropriate for detecting this channel. Refer to your instrument manual or manufacturer for support. Brilliant Violet 605™ is a trademark of Sirigen Group Ltd.</p> <p>Learn more about Brilliant Violet™.</p> <p>This product is subject to proprietary rights of Sirigen Inc. and is made and sold under license from Sirigen Inc. The purchase of this product conveys to the buyer a non-transferable right to use the purchased product for research purposes only. This product may not be resold or incorporated in any manner into another product for resale. Any use for therapeutics or diagnostics is strictly prohibited. This product is covered by U.S. Patent(s), pending patent applications and foreign equivalents.</p>
Excitation Laser	Violet Laser (405 nm)
Product Citations	<ol style="list-style-type: none">Hernández-Santos N, <i>et al.</i> 2018. Cell Host Microbe. 1.313194444. PubMedTordesillas L, <i>et al.</i> 2018. Nat Commun. 9:5238. PubMedXu W <i>et al.</i> 2019. Immunity. 50(4):1054-1068 . PubMedLong C, <i>et al.</i> 2016. Toxicol Sci. 152: 85 - 98. PubMedHuang J, <i>et al.</i> 2021. Immunity. 54(4):673-686.e4. PubMed

Antigen Details

Structure	G-protein linked, seven transmembrane, domain-spanning chemokine receptor; 367 amino acids long and about 75% homology with human CCR6
Distribution	Expressed on spleen B cells, immature splenic dendritic cells in Peyer's patches, and some T cell subset
Function	Important in regulating mucosal immune response
Ligand/Receptor	CCL20/MIP-3a
Cell Type	B cells, Dendritic cells, T cells
Biology Area	Immunology
Molecular Family	CD Molecules, Cytokine/Chemokine Receptors, GPCR
Antigen References	<ol style="list-style-type: none">1. Rossi D, <i>et al.</i> 1997. <i>J. Immunol.</i> 158:1033.2. Lukacs N, <i>et al.</i> 2001. <i>J. Exp. Med.</i> 194:551.3. Cook DN, <i>et al.</i> 2000. <i>Immunity</i> 12:495
Gene ID	12458



Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

Purified anti-mouse CD196 (CCR6), PE anti-mouse CD196 (CCR6), Alexa Fluor® 647 anti-mouse CD196 (CCR6), PerCP/Cyanine5.5 anti-mouse CD196 (CCR6), APC anti-mouse CD196 (CCR6), PE/Cyanine7 anti-mouse CD196 (CCR6), Brilliant Violet 421™ anti-mouse CD196 (CCR6), Brilliant Violet 605™ anti-mouse CD196 (CCR6), PE/Dazzle™ 594 anti-mouse CD196 (CCR6), Brilliant Violet 785™ anti-mouse CD196 (CCR6), TotalSeq™-A0225 anti-mouse CD196 (CCR6), TotalSeq™-B0225 anti-mouse CD196 (CCR6), TotalSeq™-C0225 anti-mouse CD196 (CCR6)

Product Data



C57BL/6 mouse splenocytes were stained with CD45R/B220 FITC and CCR6 (clone 29-2L17) Brilliant Violet 605™ (top) or Armenian hamster IgG Brilliant Violet 605™ isotype control (bottom).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587