

LEGENDplex™ Mouse CCL4 (MIP-1β) Capture Bead B5, 13X

Catalog# / Size	740091 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	CCL4
Description	LEGENDplex™ Mouse Proinflammatory Chemokine Mix and Match Panel Capture Beads are intended for use with the following reagents: 740074 (LEGENDplex™ Mouse Proinflammatory Chemokine Detection Antibodies) 740370 (LEGENDplex™ Mouse Proinflammatory Chemokine Standard) 740368 (LEGENDplex™ Buffer Set A)

Product Details

Verified Reactivity	Mouse
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Product Citations	1. Aguayo–Mazzucato C, <i>et al.</i> 2019. Cell Metab. 30:129. PubMed
--------------------------	---

Antigen Details

Biology Area	Cell Biology, Immunology, Innate Immunity, Neuroinflammation, Neuroscience
Molecular Family	Cytokines/Chemokines
Gene ID	20303

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

