

## Purified anti-human CD272 (BTLA) Antibody

<b>Catalog# / Size</b>	344502 / 100 µg
<b>Clone</b>	MIH26
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	BTLA, B and T lymphocyte attenuator
<b>Isotype</b>	Mouse IgG2a, κ
<b>Description</b>	B and T lymphocyte attenuator (BTLA) is an Ig superfamily coinhibitory receptor with structural similarity to programmed cell death 1 (PD-1) and CTLA-4. BTLA is expressed on B cells, T cells, macrophages, dendritic cells, NKT cells, and NK cells. Engagement of BTLA by its ligand Herpes Virus Entry Mediator (HVEM) is critical for negatively regulating immune response. The absence of BTLA with HVEM inhibitory interactions leads to increased experimental autoimmune encephalomyelitis severity, enhanced rejection of partially mismatched allografts, an increased CD8 <sup>+</sup> memory T cell population, increased severity of colitis, and reduced effectiveness of T regulatory cells. BTLA plays an important role in the induction of peripheral tolerance of both CD4 <sup>+</sup> and CD8 <sup>+</sup> T cells in vivo. Tolerant T cells have significant up-regulated expression of BTLA compared with effector and naïve T cells. BTLA may cooperate with CTLA-4 and PD-1 to control T cell tolerance and autoimmunity. It has been reported that BTLA may regulate T cell function through binding to B7-H4.

### Product Details

<b>Reactivity</b>	Human, African Green
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Human BTLA transfected cells
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	The antibody was purified by affinity chromatography.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C.
<b>Application</b>	<a href="#">FC - Quality tested</a> <a href="#">Agonist - Reported in the literature, not verified in house</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is ≤ 2.0 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.
<b>Application Notes</b>	Additional reported applications (for the relevant formats) include: inhibition of T cell proliferation and cytokine production <sup>1</sup> . Clone MIH26 has agonistic activity on BTLA, resulting in the inhibition of activation.  The Ultra-LEAF™ purified antibody (Endotoxin < 0.01 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for use in <i>in vitro</i> and <i>in vivo</i> biofunctional assays.
<b>Application References</b> (PubMed link indicates BioLegend citation)	1. Otsuki N, <i>et al.</i> 2006. <i>Biochem. Bioph. Res. Co.</i> 344:1121. 2. Okano M, <i>et al.</i> 2008. <i>Clin. Exp. Allergy</i> 38:1891.
<b>Product Citations</b>	1. Kinosada H, <i>et al.</i> 2017. <i>PLoS Pathog.</i> 10.1371/journal.ppat.1006120. <a href="#">PubMed</a> 2. Lavin Y <i>et al.</i> 2017. <i>Cell.</i> 169(4):750-765 . <a href="#">PubMed</a>
<b>RRID</b>	AB_2043950 (BioLegend Cat. No. 344502)

## Antigen Details

<b>Structure</b>	Ig superfamily, with similar structure to CTLA-4 and PD-1
<b>Distribution</b>	B cells and T lymphocytes, NK cells and dendritic cells
<b>Function</b>	Negative regulation of T cell activation, proliferation and cytokine production
<b>Ligand/Receptor</b>	HVEM
<b>Cell Type</b>	B cells, Dendritic cells, NK cells, T cells, Tregs
<b>Biology Area</b>	Immunology, Inhibitory Molecules
<b>Molecular Family</b>	CD Molecules
<b>Antigen References</b>	1. Watanabe N, <i>et al.</i> 2003. <i>Nat. Immunol.</i> 4:670. 2. Sun Y, <i>et al.</i> 2009. <i>J. Immunol.</i> 183:1946. 3. Gonzalez LC, <i>et al.</i> 2005. <i>P. Natl. Acad. Sci. USA</i> 102:1116.
<b>Gene ID</b>	<a href="#">151888</a>


## Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

## Other Formats

Purified anti-human CD272 (BTLA), PE anti-human CD272 (BTLA), APC anti-human CD272 (BTLA), Brilliant Violet 421™ anti-human CD272 (BTLA), PerCP/Cyanine5.5 anti-human CD272 (BTLA), APC/Cyanine7 anti-human CD272 (BTLA), PE/Cyanine7 anti-human CD272 (BTLA), Alexa Fluor® 647 anti-human CD272 (BTLA), PE/Dazzle™ 594 anti-human CD272 (BTLA), FITC anti-human CD272 (BTLA), TotalSeq™-A0170 anti-human CD272 (BTLA), TotalSeq™-C0170 anti-human CD272 (BTLA), Ultra-LEAF™ Purified anti-human CD272 (BTLA), PE/Cyanine5 anti-human CD272 (BTLA), TotalSeq™-B0170 anti-human CD272 (BTLA), TotalSeq™-D0170 anti-human CD272 (BTLA)

## Product Data



Human peripheral blood lymphocytes stained with CD3 APC and CD272 (BTLA) (clone MIH26) Purified (left) or Mouse IgG2a, κ Purified isotype control (right) followed by anti-mouse IgG FITC.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587