

Western-Ready™ MES SDS-PAGE Running Buffer (10X)

Catalog# / Size 426307 / 1 L
Regulatory Status RUO

Product Details

Storage & Handling Upon receipt, store solutions at room temperature. Products are shipped at ambient (room) temperature.

Application [WB - Quality tested](#)

Recommended Usage Dilute, to prepare 1X Western-Ready™ MES SDS-Page Running Buffer, add 100 mL of 10X Western-Ready™ MES SDS-Page Running Buffer to 900 mL of H₂O and mix.

Application Notes This buffer is formulated for resolving small- to medium-sized proteins on Bis-Tris SDS-PAGE gels.

Antigen Details

Gene ID NA

Related Protocols

[Western Blotting Protocol](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587