

## Western-Ready™ TBS Tween-20 Buffer (20X)

<b>Catalog# / Size</b>	426309 / 1 L
<b>Regulatory Status</b>	RUO
<b>Description</b>	This buffer is formulated to reduce non-specific binding of antibodies during immunoassay procedures such as western blotting.

### Product Details

---

<b>Formulation</b>	At 1X concentration: 137 mM NaCl, 25 mM Tris, 27 mM KCl, 0.1% Tween-20, pH 7.4
<b>Storage &amp; Handling</b>	Ambient room temperature, do not refrigerate
<b>Application</b>	<a href="#">WB - Quality tested</a>
<b>Recommended Usage</b>	To prepare 1X Western-Ready™ TBS Tween-20 Buffer, add 50 mL of Western-Ready™ TBS Tween-20 Buffer (20X) to 950 mL of H <sub>2</sub> O and mix.

### Antigen Details

---

<b>Gene ID</b>	NA
----------------	----

### Related Protocols

---

[Western Blotting Protocol](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587