

## FITC Mouse IgG3, κ Isotype Ctrl Antibody

<b>Catalog# / Size</b>	401317 / 50 µg
<b>Clone</b>	MG3-35
<b>Regulatory Status</b>	RUO
<b>Isotype</b>	Mouse IgG3, κ
<b>Description</b>	The isotype of MG3-35 immunoglobulin is mouse IgG3,κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

### Product Details

---

<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Trinitrophenol + KLH
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	The antibody was purified by affinity chromatography, and conjugated with FITC under optimal conditions.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">FC - Quality tested</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, use the isotype control at the same concentration as your primary antibody. Use our <a href="#">Concentration Lookup</a> tool to find the exact concentrations of your lots of product.
<b>Excitation Laser</b>	Blue Laser (488 nm)
<b>Application Notes</b>	The MG3-35 immunoglobulin is useful as an isotype-matched control (for the relevant formats) for Western blotting, immunoprecipitation, immunohistochemistry, functional assay, immunofluorescence microscopy, immunocytochemistry and immunofluorescent staining (surface or intracellular) for flow cytometric analysis. The LEAF™ purified antibody (Endotoxin <0.1 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for functional assays (Cat. No. 401310) as negative control.
<b>Product Citations</b>	<ol style="list-style-type: none"><li>1. Li B, <i>et al.</i> 2019. Stem Cell Res Ther. 10:257. <a href="#">PubMed</a></li><li>2. Ameres S, <i>et al.</i> 2014. J Immunol. 192:5894. <a href="#">PubMed</a></li><li>3. Bieren J, <i>et al.</i> 2015. J Immunol. 194:1154. <a href="#">PubMed</a></li><li>4. Gupta P, <i>et al.</i> 2016. PLoS One. 11: 0154725. <a href="#">PubMed</a></li></ol>
<b>RRID</b>	AB_10680461 (BioLegend Cat. No. 401317)

### Antigen Details

---

<b>Gene ID</b>	NA
----------------	----

### Other Formats

---

Purified Mouse IgG3, κ Isotype Ctrl, Biotin Mouse IgG3, κ Isotype Ctrl, FITC Mouse IgG3, κ Isotype Ctrl, Alexa Fluor® 488 Mouse IgG3, κ Isotype Ctrl, PE Mouse IgG3, κ Isotype Ctrl, Alexa Fluor® 647 Mouse IgG3, κ Isotype Ctrl, Ultra-LEAF™ Purified Mouse IgG3, κ isotype Ctrl

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587