

LEGENDplex™ Buffer Set H

Catalog# / Size	740620 / 100 tests
Regulatory Status	RUO
Description	<p>LEGENDplex™ Buffer Set H includes the following items:</p> <ul style="list-style-type: none">Setup BeadsSA-PEAssay BufferWash BufferPlate SealersMatrix AStandard Reconstitution Buffer <p>This Buffer Set does <u>not</u> include plates or the Capture Beads, Standard, and Detection Antibodies needed to complete a LEGENDplex™ Assay.</p> <p>These additional components can be purchased separately by selecting the Mix & Match tab at www.biolegend.com/legendplex, and then navigating to your desired assay.</p> <p>LEGENDplex™ Buffer Set H is designed for use with the following LEGENDplex™ Mix and Match Panels:</p> <ul style="list-style-type: none">Human Hematopoietic Stem Cell PanelHuman Vascular Inflammation Panel 2Human Proinflammatory Chemokine Panel 2

Product Details

Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Gene ID	NA
----------------	----

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

