

## Purified anti-rat CD45RA Antibody

<b>Catalog# / Size</b>	202301 / 100 µg
<b>Clone</b>	OX-33
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Leukocyte common antigen, high molecular weight form
<b>Isotype</b>	Mouse IgG1, κ
<b>Description</b>	CD45 is a protein tyrosine phosphatase with multiple isoforms differing as a result of alternative splicing of the extracellular domain and glycosylation. CD45 is expressed on all hematopoietic cells except erythrocytes and platelets. CD45RA is one of the CD45 isoforms with a molecular weight of 200-220 kD. It is expressed almost exclusively on B cells. CD45 functions in signal transduction through T and B cell antigen receptors. CD45 has been shown to interact with various proteins, including galectin-1, CD2, CD3, and CD4.

### Product Details

---

<b>Verified Reactivity</b>	Rat
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Leukocyte common antigen purified from rat splenocytes.
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	The antibody was purified by affinity chromatography.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C.
<b>Application</b>	<a href="#">FC - Quality tested</a> <a href="#">IHC-F - Reported in the literature, not verified in house</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is ≤ 0.25 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.
<b>Application Notes</b>	Additional reported applications (for the relevant formats) include: immunohistochemistry of acetone-fixed frozen sections <sup>1</sup> .
<b>Application References</b>	1. Smith T, <i>et al.</i> 2000. <i>Nature Med.</i> 6:62. (IHC)
<b>(PubMed link indicates BioLegend citation)</b>	
<b>RRID</b>	AB_314007 (BioLegend Cat. No. 202301)

### Antigen Details

---

<b>Structure</b>	Protein tyrosine phosphatase, type I transmembrane protein, high molecular weight isoform approximately 200-220 kD
<b>Distribution</b>	Primarily on B lymphocytes
<b>Function</b>	B cell receptor signal transduction
<b>Ligand/Receptor</b>	Galectin-1, CD2, CD3, CD4

<b>Cell Type</b>	B cells
<b>Biology Area</b>	Cell Biology, Immunology, Neuroscience, Neuroscience Cell Markers, Inhibitory Molecules
<b>Molecular Family</b>	CD Molecules
<b>Antigen References</b>	1. Sunderland CA, <i>et al.</i> 1979. <i>Eur. J. Immunol.</i> 9:155. 2. Woolett GR, <i>et al.</i> 1985. <i>Eur. J. Immunol.</i> 15:168.
<b>Gene ID</b>	<a href="#">24699</a>


## Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

## Other Formats

Purified anti-rat CD45RA, Biotin anti-rat CD45RA, FITC anti-rat CD45RA, PE anti-rat CD45RA, Alexa Fluor® 488 anti-rat CD45RA, APC anti-rat CD45RA, PerCP/Cyanine5.5 anti-rat CD45RA, PE/Cyanine7 anti-rat CD45RA

## Product Data



LOU rat splenocytes stained with purified OX-33, followed by anti-mouse IgG FITC

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587