

Recombinant Mouse IL-3 (carrier-free)

Catalog# / Size	575502 / 10 µg 575504 / 25 µg 575506 / 100 µg 575508 / 500 µg
Regulatory Status	RUO
Other Names	Interleukin-3, Burst promoting activity, Eosinophil colony stimulating factor (Eo-CSF), Hematopoietic cell growth factor (HCGF), Mast (MGF/MCGF), Multi-colony stimulating (Multi-CSF), P cell stimulating activity (PCSA), Thy1 inducing factor
Description	IL-3 is the most potent growth factor for basophils followed by granulocyte-macrophage colony-stimulating factor and IL-5. These cytokines also act on mature basophils through specific receptors, thereby mediating adhesion, migration, and releasability. IL-3 is highly expressed by mast cells, and rapid and large amount of autocrine IL-3 production is responsible for mast cell survival by IgE in the absence of antigen. IL-3 has also been implicated in the pathogenesis of several chronic inflammatory diseases, including asthma, atherosclerosis, and neurodegenerative disorders, such as multiple sclerosis.


Product Details

Source	Mouse IL-3, amino acids Ala 27 - Cys 166 (Accession# NM_010556), was expressed in <i>E. coli</i> .
Molecular Mass	The 140 amino acid recombinant protein has a predicted molecular mass of 15673.7 Da. The DTT-reduced protein migrates at approximately 16 kDa and the non-reduced protein migrates with slightly greater mobility by SDS-PAGE.
Purity	>98%, as determined by Coomassie stained SDS-PAGE.
Formulation	0.22 µm filtered protein solution is in 10 mM NaHPO ₄ pH 7.2, 0.15 M NaCl.
Endotoxin Level	Less than 0.01ng per µg cytokine as determined by the LAL method.
Concentration	10 and 25 µg sizes are bottled at 200 µg/mL. 100 µg size and larger sizes are lot-specific and bottled at the concentration indicated on the vial. To obtain lot-specific concentration, please enter the lot number in our Concentration and Expiration Lookup or Certificate of Analysis online tools.
Storage & Handling	Unopened vial can be stored between 2°C and 8°C for up to 2 weeks, at -20°C for up to six months, or at -70°C or colder until the expiration date. For maximum results, quick spin vial prior to opening. The protein can be aliquoted and stored at -20°C or colder. Stock solutions can also be prepared at 50 - 100 µg/mL in appropriate sterile buffer, carrier protein such as 0.2 - 1% BSA or HSA can be added when preparing the stock solution. Aliquots can be stored between 2°C and 8°C for up to one week and stored at -20°C or colder for up to 3 months. Avoid repeated freeze/thaw cycles.
Activity	ED ₅₀ = 0.05 - 0.200 ng/ml, corresponding to a specific activity of 0.5-2.0 x 10 ⁷ units/mg, as determined by the dose dependent stimulation of a M-NFS-60 cell proliferation assay.
Application	Bioassay
Application Notes	BioLegend carrier-free recombinant proteins provided in liquid format are shipped on blue-ice. Our comparison testing data indicates that when handled and stored as recommended, the liquid format has equal or better stability and shelf-life compared to commercially available lyophilized proteins after reconstitution. Our liquid proteins are verified in-house to maintain activity after shipping on blue ice and are backed by our 100% satisfaction guarantee . If you have any concerns, contact us at tech@biolegend.com .
Additional Product Notes	View more applications data for this product in our Scientific Poster Library .
Product Citations	1. Paris J <i>et al.</i> 2019. Cell Stem Cell. 25(1):137-148 . PubMed 2. Sinclair JKL, <i>et al.</i> 2018. Cell Chem Biol. 25:857. PubMed

Antigen Details

Structure	Cytokine
Distribution	Activated T cells, Th1, Th2, mast cells, eosinophils, keratinocytes, NK cells, and endothelial cells.
Function	Interleukine-3 is a pleiotropic factor which promotes both self-renewal and differentiation of early multipotential progenitors and is involved in inducible hematopoiesis in response to infections. IL-3 stimulates colony formation of megakaryocytes, neutrophils, and macrophages from bone marrow cultures.
Interaction	Erythroid cells, megakaryocytes, neutrophils, eosinophils, basophils, mast cells, monocytic lineages
Ligand/Receptor	Heterodimer IL-3R α (CD123) and two alternate β -subunits, AIC2A binds IL-3 and is specific to IL-3R α , AIC2B (CDw131) shared in common with IL-5R, GM-CSFR.
Cell Type	Hematopoietic stem and progenitors
Biology Area	Stem Cells
Molecular Family	Cytokines/Chemokines
Antigen References	<ol style="list-style-type: none">1. Kohno M, <i>et al.</i> 2005 <i>Blood</i> 105:2059-2065.2. Valent P and Dahinden CA, <i>et al.</i> 2010 <i>Curr Opin Hematol</i> 17:60-66.3. Kleemann R, <i>et al.</i> 2008 <i>Cardiovascular Research</i> 79:360-376.4. Murphy JM and Young IG, <i>et al.</i> 2006 <i>Vitam Horm</i> 74:1-30.
Gene ID	16187

Product Data



M-NFS-60 cell proliferation induced by mouse IL-3.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587