

LEGENDplex™ Human TNF-α Capture Beads A5, 13X

Catalog# / Size	740513 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	TNF-alpha TNF-a, TNF-alpha, TNFa, TNFalpha, TNFSF2, TNLG1F, Tumor necrosis factor alpha, TNF alpha, cachexin, cachectin
Description	<p>The LEGENDplex™ Human Macrophage/Microglia Panel Mix and Match Beads are individual capture beads for 13-targets including IL-12p70, TNF-α, IL-6, IL-4, IL-10, IL-1β, Arginase, TARC, IL-1RA, IL-12p40, IL-23, IFN-γ, IP-10. The beads are recommended for use with each other and the following in Mix and Match Assays:</p> <p>740511 (LEGENDplex™ Human Macrophage/Microglia Panel Standard) 740510 (LEGENDplex™ Human Macrophage/Microglia Panel Detection Antibodies) 741149 (LEGENDplex™ Buffer Set A6)</p>

Product Details

Verified Reactivity	Human
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Cell Biology, Immunology, Neuroinflammation, Neuroscience
Molecular Family	Cytokines/Chemokines
Gene ID	7124

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

