

Purified anti-Histone H3 (C-terminus) Antibody

Catalog# / Size	819411 / 25 µg 819412 / 100 µg
Clone	1B1-B2
Regulatory Status	RUO
Other Names	Histone H3.1, histone H3/a, histone 1, H3a, H3 histone family, member A, H3/A, H3FA
Isotype	Mouse IgG3, κ
Description	<p>Histone proteins are classified into core histones (H2A, H2B, H3, H4) and linker histones (H1, H5). Core histones form an octamer, which contains two H2A-H2B dimers and one H3-H4 tetramer. Core histones are predominantly globular except for the unstructured N-terminal tails. Posttranslational modifications, such as acetylation, methylation, phosphorylation, ubiquitination, SUMOylation and ADP-ribosylation occur in histone tails.</p> <p>Histone modifications induce changes of chromatin structure and thereby affect the accessibility of transcription factors, nuclear proteins and enzymes to genomic DNA, resulting in gene activation or repression. It is known that histone modifications play critical roles in DNA repair, DNA replication, transcription regulation, alternative splicing and chromosome condensation and some diseases including autoimmune diseases and cancers.</p>

Product Details

Verified Reactivity	Human, Mouse, Rat, Drosophila
Antibody Type	Monoclonal
Host Species	Mouse
Immunogen	This monoclonal antibody was raised against a peptide sequence derived from the C-terminal domain of human Histone H3 protein.
Formulation	50 mM sodium citrate, 150mM sodium chloride, pH 5.8 with stabilizer
Preparation	The antibody was purified by affinity chromatography.
Concentration	0.5 mg/mL
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C.
Application	IHC-P - Quality tested WB - Verified
Recommended Usage	Each lot of this antibody is quality control tested by immunohistochemistry. For immunohistochemistry, a concentration of 1.0 µg/mL is suggested. For Western blotting, the suggested use of this reagent is 0.01 - 0.5 µg/mL. It is recommended that the reagent be titrated for optimal performance for each application.
Application Notes	Reactivity to Drosophila was only verified with the purified format.
Additional Product Notes	For Western blotting, the suggested use of this reagent is 0.01 - 0.2 µg/mL in human, mouse, and rat and 0.1 - 0.5 µg/mL in Drosophila.
Product Citations	1. Bado IL, <i>et al.</i> 2021. Developmental Cell. 56(8):1100-1117.e9. PubMed
RRID	AB_2820127 (BioLegend Cat. No. 819411) AB_2820128 (BioLegend Cat. No. 819412)

Antigen Details

Structure	Histone proteins H3 and H4 bind to form a H3-H4 dimer, two of these H3-H4 dimers combine to form a tetramer. This tetramer further combines with two H2a-H2b dimers to form the compact Histone octamer core.
Distribution	Nucleus
Function	Histones play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.
Interaction	Rsk-2, MSK1, GCN5 family of HATs, HDACs, PARP, SET7/9, and CARM1
Biology Area	Cell Biology, Chromatin Remodeling/Epigenetics, Neuroscience, Transcription Factors
Antigen References	1. Chen HM, <i>et al.</i> 2015. <i>Epigenetics and Dermatology</i> 409.
Gene ID	8350

Related Protocols


[Western Blotting Protocol](#)

[Immunohistochemistry Protocol for Paraffin-Embedded Sections](#)


Other Formats

Alexa Fluor® 594 anti-Histone H3 (C-terminus), HRP anti-Histone H3 (C-terminus), Alexa Fluor® 647 anti-Histone H3 (C-terminus), Purified anti-Histone H3 (C-terminus), Go-ChIP-Grade™ Purified anti-Histone H3 (C-terminus)


Product Data




IHC staining of purified anti-Histone H3 (C-Terminus) antibody (clone 1B1-B2) on formalin-fixed paraffin-embedded rat brain tissue. Following antigen retrieval using Sodium Citrate H.I.E.R., the tissue was incubated with 1 µg/mL of the primary antibody overnight at 4°C. BioLegend's Ultra Streptavidin (USA) HRP kit (Multi-Species, DAB, Cat. No. 929901) was used for detection followed by hematoxylin counterstaining, according to the protocol provided. The image was captured with a 40X objective.




IHC staining of purified anti-Histone H3 (C-Terminus) antibody (clone 1B1-B2) on formalin-fixed paraffin-embedded mouse brain tissue. Following antigen retrieval using Sodium Citrate H.I.E.R., the tissue was incubated with 1 µg/mL of the primary antibody overnight at 4°C. BioLegend's Ultra Streptavidin (USA) HRP kit (Multi-Species, DAB, Cat. No. 929901) was used for detection followed by hematoxylin counterstaining, according to the protocol provided. The image was captured with a 40X objective. Scale bar: 50 µm



IHC staining of purified anti-Histone H3 (C-Terminus) antibody (clone 1B1-B2) on formalin-fixed paraffin-embedded rat brain tissue. Following antigen retrieval using Sodium Citrate H.I.E.R., the tissue was incubated with 1 µg/mL of the primary antibody overnight at 4°C. BioLegend's Ultra Streptavidin (USA) HRP kit (Multi-Species, DAB, Cat. No. 929901) was used for detection followed by hematoxylin counterstaining, according to the protocol provided. The image was captured with a 40X objective. Scale bar: 50 µm



Western blot of purified anti-Histone H3 (C-terminus) antibody (clone 1B1-B2). Lane 1: Molecular weight marker; Lane 2: 20 µg of human brain lysate; Lane 3: 20 µg of mouse brain lysate; Lane 4: 20 µg of rat brain lysate; Lane 5: 10 µg of Drosophila head lysate. The blot was incubated with 0.1 µg/mL of the primary antibody overnight at 4°C, followed by incubation with HRP Goat anti-mouse IgG Antibody (Cat. No. 405306). Enhanced chemiluminescence was used as the detection system.



Western blot of purified anti-Histone H3 (C-terminus) antibody (clone 1B1-B2). Lane 1: Molecular weight marker; Lane 2: 20 µg of Drosophila head lysate; Lane 3: 20 µg of Drosophila S2 (embryonic) cell lysate. The blot was incubated with 0.1 µg/mL of the primary antibody overnight at 4°C, followed by incubation with HRP Goat anti-mouse IgG Antibody (Cat. No. 405306). Enhanced chemiluminescence was used as the detection system.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587