

Brilliant Violet 421™ anti-human CD200R Antibody

Catalog# / Size	329313 / 25 tests 329314 / 100 tests
Clone	OX-108 (See other available formats)
Other Names	Orexin Receptor 2 (OX2R), CD200 Receptor
Isotype	Mouse IgG1, κ
Description	CD200R, also known as OX2R, is a membrane glycoprotein with up to 70% of its weight derived from N-linked glycosylation. CD200R is expressed primarily by monocytes and neutrophils, but also by other leukocytes including T cells and mast cells. The interaction between CD200 and CD200R may contribute to pathways that suppress and limit macrophage induced inflammatory damage in tissues.

Product Details

Reactivity	Human, Baboon
Antibody Type	Monoclonal
Host Species	Mouse
Immunogen	Human CD200R full length fusion protein
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA).
Preparation	The antibody was purified by affinity chromatography and conjugated with Brilliant Violet 421™ under optimal conditions. The solution is free of unconjugated Brilliant Violet 421™ and unconjugated antibody.
Concentration	Lot-specific (please contact technical support for concentration and total µg amount, or use our Lookup tool if you have a lot number.)
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	<p>Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is 5 µl per million cells in 100 µl staining volume or 5 µl per 100 µl of whole blood.</p> <p>Brilliant Violet 421™ excites at 405 nm and emits at 421 nm. The standard bandpass filter 450/50 nm is recommended for detection. Brilliant Violet 421™ is a trademark of Sirigen Group Ltd.</p> <p>Learn more about Brilliant Violet™.</p> <p>This product is subject to proprietary rights of Sirigen Inc. and is made and sold under license from Sirigen Inc. The purchase of this product conveys to the buyer a non-transferable right to use the purchased product for research purposes only. This product may not be resold or incorporated in any manner into another product for resale. Any use for therapeutics or diagnostics is strictly prohibited. This product is covered by U.S. Patent(s), pending patent applications and foreign equivalents.</p>
Excitation Laser	Violet Laser (405 nm)
Application References	1. Wright GJ, et al. 2003. J. Immunol. 171:3034.
(PubMed link indicates BioLegend citation)	

Antigen Details

Structure	Immunoglobulin superfamily glycoprotein.
Distribution	Monocytes, neutrophils, T cells, mast cells.
Ligand Receptor	CD200, OX2.

Cell Type	Mast cells, Monocytes, Neutrophils, T cells
Biology Area	Cell Biology, Immunology, Innate Immunity, Neuroscience, Neuroscience Cell Markers
Molecular Family	CD Molecules
Antigen References	1. Wright GJ, et al. 2003. J. Immunol. 171:3034. 2. Gorczynski R, et al. 2004. J. Immunol. 172:7744.
Gene ID	131450


Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

Purified anti-human CD200R, PE anti-human CD200R, APC anti-human CD200R, PE/Dazzle™ 594 anti-human CD200R, PE/Cy7 anti-human CD200R

Product Data



Human peripheral blood lymphocytes were stained with CD3 (clone UCHT1) PE and CD200R (clone OX-108) Brilliant Violet 421™ (left) or Mouse IgG1, κ Brilliant Violet 421™ isotype control (right).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 9727 Pacific Heights Blvd, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587