

## LEGENDplex™ Human EGF Capture Bead A7, 13X

<b>Catalog# / Size</b>	740703 / 100 tests (270 µL)
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Epidermal Growth Factor, Epidermal Growth Factor (Beta-Urogastrone), Pro-Epidermal Growth Factor, Beta-Urogastrone, HOMG4, URG
<b>Description</b>	<p>The LegendPlex™ Human Angiogenesis Panel 1 V02 Mix and Match Subpanel includes beads for 10 analytes: IL-6, Angiopoietin-1, Angiopoietin-2, EGF, FGF-basic, CXCL8 (IL-8), PECAM-1, PIGF, VEGF, TNF-<math>\alpha</math>. These beads are recommended for use with each other and the following reagent in a Mix and Match Assay:</p> <p>741216 LEGENDplex™ Human Angiogenesis Panel 1 Detection Abs V02 740699 LEGENDplex™ Human Angiogenesis Panel 1 Standard V02 741218 LEGENDplex™ Buffer Set I1 740377 LEGENDplex™ Filter Plate OR 740379 LEGENDplex™ V-bottom Plate</p>

### Product Details

---

<b>Verified Reactivity</b>	Human
<b>Application</b>	<a href="#">Multiplex</a>
	Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a>
	Download the <a href="#">LEGENDplex™ software here</a> .

### Antigen Details

---

<b>Biology Area</b>	Cell Biology, Immunology, Innate Immunity, Neuroscience, Stem Cells, Synaptic Biology
<b>Molecular Family</b>	Cytokines/Chemokines, Growth Factors
<b>Gene ID</b>	<a href="#">1950</a>

### Related Protocols


---

[LEGENDplex™ Laboratory Protocol - Video](#)

### Product Data

---


**LEGENDplex™ Human EGF Capture Bead A7  
Representative Standard Curve**



- Data generated using the following LEGENDplex™ Mix & Match Reagents:
- Human EGF Capture Bead A7, 13X (Cat. # 740703)
  - Human Angiogenesis Panel 1 Detection Antibodies (Cat. # 740700)
  - Human Angiogenesis Panel 1 Standard (Cat. # 740699)
  - Buffer Set I (Cat. # 740623)

**Note: Data acquired on a BD FACSCalibur™ flow cytometer.**

**EGF**



Human plasma samples (n = 80) were quantified using the following LEGENDplex™ Mix & Match Reagents:

- Human EGF Capture Bead A7, 13X (Cat. # 740703)
- Human Angiogenesis Panel 1 Detection Antibodies (Cat. # 740700)
- Human Angiogenesis Panel 1 Standard (Cat. # 740699)
- Buffer Set I (Cat. # 740623)

Data are presented as Mean + SD. Plasma samples were provided courtesy of Professor Huaizhu Wu, Baylor College of Medicine.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587