

LEGENDplex™ Human IL-10 Capture Bead B7, 13X

Catalog# / Size	740358 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	Interleukin-10 (IL10), CSIF, B-TCGF, BTCGF, TGIF
Description	The Human IL-10 Capture Bead is suitable for use with the following Mix and Match panels:

Human Anti-Virus Response Panel

740364 (LEGENDplex™ Human Anti-Virus Response Panel Detection Abs)
740365 (LEGENDplex™ Human Anti-Virus Response Panel Standard)
740368 (LEGENDplex™ Buffer Set A)

Human Vascular Inflammation Panel 2

740968 (LEGENDplex™ Human Vascular Inflammation Panel 2 Detection Abs)
740967 (LEGENDplex™ Human Vascular Inflammation Panel 2 Standard)
740620 (LEGENDplex™ Buffer Set H)

For product manuals for these panels, please visit the Mix and Match tab on the [LEGENDplex™](#) page.

Product Details

Verified Reactivity	Human
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Immunology
Molecular Family	Cytokines/Chemokines
Gene ID	NA

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587