

LEGENDplex™ Human IL-4 Capture Beads A6, 13X

Catalog# / Size	740540 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	IL-4, IL4, BSF1, IL-4, Interleukin 4, BCGF1, BSF-1, BCGF-1
Description	LEGENDplex™ Human B Cell Panel Capture Beads are intended for use with the following reagents: 740536 (LEGENDplex™ Human B Cell Panel Detection Antibodies) 740537 (LEGENDplex™ Human B Cell Panel Standard) 741150 (LEGENDplex™ Buffer Set A7)

Product Details

Verified Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Immunology
Molecular Family	Cytokines/Chemokines
Gene ID	3565

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587