

LEGENDplex™ Mouse IL-1β Capture Bead B2, 13X

Catalog# / Size	740157 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	IL-1beta, IL-1
Description	This Mouse IL-1β Capture Bead is suitable for use with the following Mix and Match panels:

Mouse Inflammation Panel

740165 (LEGENDplex™ Mouse Inflammation Panel Detection Antibodies)
740371 (LEGENDplex™ Mouse Inflammation Panel Standard)
740373 (LEGENDplex™ Buffer Set B)

Mouse Anti-Virus Response Panel

740625 (LEGENDplex™ Mouse Anti-Virus Response Panel Detection Antibodies)
740624 (LEGENDplex™ Mouse Anti-Virus Response Panel Standard)
740623 (LEGENDplex™ Buffer Set I)

For product manuals for these panels, please visit the Mix and Match tab on the [LEGENDplex™](#) page.

Product Details

Verified Reactivity	Mouse
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .
Product Citations	1. Meek SM, <i>et al.</i> 2018. Pathogens. 7:22. PubMed

Antigen Details

Biology Area	Cell Biology, Immunology, Innate Immunity, Neuroinflammation, Neuroscience
Molecular Family	Cytokines/Chemokines
Gene ID	16176

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587