

LEGENDplex™ Rat IL-6 Capture Bead B9, 13X

Catalog# / Size	740264 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	Interleukin-6, Interferon-beta2, B cell stimulating factor-2 (BSF-2), Cytotoxic T cell differentiation factor (CDF), Hepatocyte stimulating factor (HSF), Hybridoma/plasmacytoma growth factor (HPGF)
Description	<p>LEGENDplex™ Rat Inflammation Mix and Match Panel Capture Beads are intended for use with the following reagents:</p> <p>740266 (LEGENDplex™ Rat Inflammation Panel Detection Antibodies) 740265 (LEGENDplex™ Rat Inflammation Panel Standard) 740373 (LEGENDplex™ Buffer Set B)</p> <p>For product manuals for these panels, please visit the Mix and Match tab on the LEGENDplex™ page.</p>

Product Details

Verified Reactivity	Rat
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Immunology, Innate Immunity, Stem Cells
Molecular Family	Cytokines/Chemokines
Gene ID	24498

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

