

## Alexa Fluor® 488 anti-Bcl-2 Antibody

<b>Catalog# / Size</b>	658703 / 25 tests 658704 / 100 tests
<b>Clone</b>	100
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	B-Cell CLL/Lymphoma 2 (Bcl-2), Protein Phosphatase 1, Regulatory Subunit 50
<b>Isotype</b>	Mouse IgG1
<b>Description</b>	Bcl-2 is a conserved anti-apoptotic protein that plays important roles in normal immunity. It forms homodimers or heterodimers with other Bcl-2 family members. This protein regulates apoptosis through controlling mitochondrial fusion and fission. Phosphorylation of Bcl-2 has been shown to enhance activity to allow response to extracellular growth-factor-mediated signals. Bcl-2 is also regulated by IRES, miR-15a, miR-16-1 and RNA-BP nucleolin. Chromosome 18q21.3 translocation and overexpression of Bcl-2 is frequently observed in follicular lymphomas and some diffuse large B-cell lymphomas.

### Product Details

---

<b>Reactivity</b>	Human
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Synthetic peptide of amino acids 41-54 of human bcl-2 oncoprotein
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)
<b>Preparation</b>	The antibody was purified by affinity chromatography and conjugated with Alexa Fluor® 488 under optimal conditions.
<b>Concentration</b>	Lot-specific (to obtain lot-specific concentration, please enter the lot number in our <a href="#">Concentration and Expiration Lookup</a> or <a href="#">Certificate of Analysis</a> online tools.)
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">ICFC - Quality tested</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">intracellular immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is 5 µl per million cells in 100 µl staining volume or 5 µl per 100 µl of whole blood.  * Alexa Fluor® 488 has a maximum emission of 519 nm when it is excited at 488 nm.  Alexa Fluor® and Pacific Blue™ are trademarks of Life Technologies Corporation.  <a href="#">View full statement regarding label licenses</a>
<b>Excitation Laser</b>	Blue Laser (488 nm)
<b>Application References</b>	<ol style="list-style-type: none"><li>1. Pezzella F, <i>et al.</i> 1993. <i>New Eng. J. Med.</i> 329:690. (IHC)</li><li>2. Pezzella F, <i>et al.</i> 1990. <i>Am. J. Pathol.</i> 137:225. (IHC)</li><li>3. Radtke AJ, <i>et al.</i> 2020. <i>Proc Natl Acad Sci U S A.</i> 117:33455-65. (SB) <a href="#">PubMed</a></li><li>4. Radtke AJ, <i>et al.</i> 2022. <i>Nat Protoc.</i> 17:378-401. (SB) <a href="#">PubMed</a></li></ol>
<b>(PubMed link indicates BioLegend citation)</b>	
<b>Product Citations</b>	<ol style="list-style-type: none"><li>1. Ramos CV, <i>et al.</i> 2020. <i>Cell Reports.</i> 32(3):107910. <a href="#">PubMed</a></li></ol>
<b>RRID</b>	AB_2563151 (BioLegend Cat. No. 658703) AB_2563152 (BioLegend Cat. No. 658704)

## Antigen Details

---

<b>Structure</b>	Two isoforms, 239 or 205 amino acids in length, with a predicted molecular weight of 26 kD or 22 kD and four 'BH' domains.
<b>Distribution</b>	Outer mitochondrial membrane, nuclear membrane, endoplasmic reticulum membrane.
<b>Function</b>	Regulates apoptosis through controlling mitochondrial fusion and fission.
<b>Interaction</b>	BAX, BAD, BAK, Bcl-X(L), EI24, APAF1, BBC3, BCL2L1, BNIP1, MRPL41, TP53BP2, FKBP8 BAG1, RAF1, EGLN3, G0S2, and BOP.
<b>Cell Type</b>	B cells
<b>Biology Area</b>	Apoptosis/Tumor Suppressors/Cell Death, Cell Biology, Cell Cycle/DNA Replication, Chromatin Remodeling/Epigenetics, Immunology, Mitochondrial Function, Neuroscience, Signal Transduction, Transcription Factors, Ubiquitin/Protein Degradation
<b>Antigen References</b>	<ol style="list-style-type: none"><li>1. Lin P, <i>et al.</i> 2013. <i>Curr. Hematol. Malig. Rep.</i> 8:243.</li><li>2. Tomita N. 2011. <i>J. Clin. Exp. Hematop.</i> 51:7.</li><li>3. Kelly PN, <i>et al.</i> 2011. <i>Cell Death Differ.</i> 18:1414.</li><li>4. Willmott S, <i>et al.</i> 2010. <i>Biochem. Soc. Trans.</i> 38:1571.</li><li>5. Rolland SG, <i>et al.</i> 2010. <i>Curr. Opin. Cell Biol.</i> 22:852.</li></ol>
<b>Gene ID</b>	<a href="#">596</a>

## Related Protocols

---

[Intracellular Cytokine Staining Protocol - Video](#)

[Intracellular Flow Cytometry Staining Protocol](#)


## Other Formats

---

Purified anti-Bcl-2, Alexa Fluor® 647 anti-Bcl-2, Alexa Fluor® 488 anti-Bcl-2, PE anti-Bcl-2, Brilliant Violet 421™ anti-Bcl-2

## Product Data

---



Human peripheral blood lymphocytes were treated with Fixation Buffer (Cat# 420801) and Permeabilization Wash Buffer (Cat# 421002), then stained with Bcl-2 (clone 100) Alexa Fluor® 488 (filled histogram) or mouse IgG1, κ Alexa Fluor® 488 isotype control (open histogram).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587

