

## Purified anti-STAT3 Antibody

<b>Catalog# / Size</b>	691502 / 100 µg
<b>Clone</b>	P83B5
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Signal transducer and activator of transcription 3, Acute-phase response factor (APRF), HIES, ADMIO
<b>Isotype</b>	Mouse IgG1, κ
<b>Description</b>	STAT3 is an 88 kD member of the STAT (signal transducer and activators of transcription) protein family that is phosphorylated in response to a cytokine receptor-associated kinase activity. Phosphorylation of STAT3 induces nuclear translocation to activate transcription. STAT3 forms both homo- and heterotrimers and is involved in the activation of genes required for cell growth and apoptosis. STAT3 is also involved in gp130 signaling and binds to IL-6 response elements in various acute phase protein promoters. STAT3 is phosphorylated by signaling from IFNs, EGF, FGF, IL-5, HGF, LIF, and BMP2. STAT3 activity is inhibited by PIAS3 and GRIM-19 and can also be regulated by the Rac1 protein.

### Product Details

---

<b>Verified Reactivity</b>	Human, Mouse
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Purified recombinant fragment (aa 20-240) of human STAT3 expressed in <i>E. coli</i> .
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	The antibody was purified by affinity chromatography.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C.
<b>Application</b>	<a href="#">WB - Quality tested</a> <a href="#">ChIP - Verified</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">Western blotting</a> . For Western blotting, the suggested use of this reagent is 0.5 - 2.0 µg per ml. For ChIP applications, the suggested use of this reagent is 1-50-1:150 dilution by volume. It is recommended that the reagent be titrated for optimal performance for each application.
<b>Application Notes</b>	The canonical sequence of STAT3 contains 770 amino acids with a predicted molecular weight of 88 kD. A shorter isoform of STAT3 with 722 amino acids (83 kD) was reported. Clone P83B5 may recognize both isoforms.
<b>RRID</b>	AB_2632650 (BioLegend Cat. No. 691502)

### Antigen Details

---

<b>Structure</b>	770 amino acids with a predicted molecular weight of 88.1 kD that contains a C-terminal SH2 domain responsible for homo- or heterodimerization.
<b>Distribution</b>	Cytoplasm.
<b>Function</b>	Phosphorylated in response to cytokine signaling by receptor-associated kinases, translocates to the nucleus to act as a transcriptional activator, activates a wide variety of genes involved in cell growth and apoptosis, and is involved in gp130-related signaling.

<b>Interaction</b>	STAT3 forms a homodimer or a heterodimer with other STAT family members and interacts with IL31RA, NCOA1, PELP1, SIPAR, SOCS7, STATIP1, CAV2, DAPK3, EIF2AK2, PIAS3, and TMF1.
<b>Cell Type</b>	Embryonic Stem Cells, Mesenchymal Stem Cells
<b>Biology Area</b>	Cell Biology, Neuroscience, Neuroscience Cell Markers, Stem Cells, Transcription Factors
<b>Molecular Family</b>	Nuclear Markers
<b>Antigen References</b>	<ol style="list-style-type: none"> <li>1. Lufe C, <i>et al.</i> 2003. <i>EMBO J.</i> 22:1325.</li> <li>2. Deo D, <i>et al.</i> 2002. <i>J. Biol. Chem.</i> 277:21237.</li> <li>3. Pfeffer L, <i>et al.</i> 1997. <i>Science.</i> 276:1418.</li> <li>4. Akira S, <i>et al.</i> 1994. <i>Cell.</i> 77:63.</li> </ol>
<b>Gene ID</b>	<a href="#">6774</a>

## Related Protocols

[BioLegend's Tools for Chromatin Immunoprecipitation \(ChIP\) Assays - Video](#)


[Chromatin Immunoprecipitation \(ChIP\) Assay Protocol](#)

[Western Blotting Protocol](#)


## Other Formats

Purified anti-STAT3, Go-ChIP-Grade™ Purified anti-STAT3

## Product Data



Total cell lysates (15 µg protein) from Raw264.7 (lane 1), PC3 (lane 2) and HeLa (lane 3) and Molt-4 (lane 4) cells were resolved by electrophoresis (4-12% Bis-Tris gel), transferred to nitrocellulose, and probed with 1 µg/mL (1:500 dilution) of Purified anti-STAT3 Antibody (clone P83B5) (upper panel). Proteins were visualized by chemiluminescence detection using a 1:3000 diluted goat anti-mouse-IgG secondary antibody conjugated to HRP for anti-STAT3 Antibody or 1:5000 diluted Direct-Blot™ HRP anti-β-Actin antibody, clone 2F1-1(lower panel) as a loading control. Lane M: Molecular weight ladder.



Chromatin Immunoprecipitation (ChIP) was performed using commercial Protein-G coated 96 well high-throughput ChIP assay kit by loading 3 µg of cross-linked chromatin samples from HeLa cells starved overnight and then treated with IL-6 with either A) 1:500 dilution of Go-ChIP-Grade™ purified anti-STAT3 (clone P83B5), or B) equal amount of purified Mouse IgG1, κ isotype control antibody. The enriched DNA was purified and quantified by real-time qPCR using primers targeting human IRF1 gene region. The amount of immunoprecipitated DNA in each sample is represented as signal relative to the 5% of total amount of input chromatin.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587