

APC/Fire™ 750 anti-mouse CD279 (PD-1) Antibody

Catalog# / Size	135239 / 25 µg 135240 / 100 µg
Clone	29F.1A12
Regulatory Status	RUO
Other Names	PD-1, Programmed Death-1
Isotype	Rat IgG2a, κ
Description	CD279, also known as programmed death-1 (PD-1), is a 50-55 kD glycoprotein belonging to the CD28 family of the Ig superfamily. PD-1 is expressed on activated splenic T and B cells and thymocytes. It is induced on activated myeloid cells as well. PD-1 is involved in lymphocyte clonal selection and peripheral tolerance through binding its ligands, B7-H1 (PD-L1) and B7-DC (PD-L2). It has been reported that PD-1 and PD-L1 interactions are critical to positive selection and play a role in shaping the T cell repertoire. PD-L1 negative costimulation is essential for prolonged survival of intratesticular islet allografts.

Product Details

Verified Reactivity	Mouse
Antibody Type	Monoclonal
Host Species	Rat
Immunogen	PD-1 cDNA followed by PD-1-Ig fusion protein
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Preparation	The antibody was purified by affinity chromatography and conjugated with APC/Fire™ 750 under optimal conditions.
Concentration	0.2 mg/ml
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis . For flow cytometric staining, the suggested use of this reagent is ≤0.25 µg per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application. * APC/Fire™ 750 has a maximum excitation of 650 nm and a maximum emission of 787 nm.
Application Notes	Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen tissue ³ , <i>in vivo</i> blocking of PD-1 binding to its ligands ^{2,3} , and spatial biology (IBEX) ^{5,6} .
Application References	<ol style="list-style-type: none">1. Good-Jacobson KL, <i>et al.</i> 2010. <i>Nat. Immunol.</i> 11:535. (FC) PubMed2. Lázár-Molnár E, <i>et al.</i> 2008. <i>Proc. Natl. Acad. Sci. USA</i> 105:2658. (Block)3. Liang SC, <i>et al.</i> 2003. <i>Eur. J. Immunol.</i> 33:2706. (FC, IHC, Block)4. Tobias J, <i>et al.</i> 2020. <i>Front Immunol.</i> 11:895 (FC, ELISA) PubMed5. Radtke AJ, <i>et al.</i> 2020. <i>Proc Natl Acad Sci U S A.</i> 117:33455-65. (SB) PubMed6. Radtke AJ, <i>et al.</i> 2022. <i>Nat Protoc.</i> 17:378-401. (SB) PubMed
Product Citations	<ol style="list-style-type: none">1. Tomaru U <i>et al.</i> 2019. <i>Cell reports.</i> 26(3):639-651 . PubMed2. Levesque S, <i>et al.</i> 2019. <i>Oncoimmunology.</i> 8:e1657375. PubMed3. Buqué A, <i>et al.</i> 2020. <i>Nat Commun.</i> 3.110416667. PubMed4. Abdelsamed HA, <i>et al.</i> 2020. <i>Nat Immunol.</i> 1.276388889. PubMed

RRID AB_2629767 (BioLegend Cat. No. 135239)
AB_2629768 (BioLegend Cat. No. 135240)

Antigen Details

Structure	A 50-55 kD glycoprotein belonging to the CD28 family of the Ig superfamily.
Distribution	Induced on splenic T and B lymphocytes, thymocytes, and myeloid cells after stimulation.
Function	Involved in lymphocyte clonal selection and peripheral tolerance, prolonged survival of allografts.
Ligand/Receptor	B7-H1 (PD-L1) and B7-DC (PD-L2)
Cell Type	B cells, T cells
Biology Area	Cancer Biomarkers, Immunology, Inhibitory Molecules
Molecular Family	CD Molecules, Immune Checkpoint Receptors
Antigen References	<ol style="list-style-type: none">1. Nishimura H, <i>et al.</i> 2001. <i>Science</i> 291:3192. Agata Y, <i>et al.</i> 1996. <i>Int. Immunol.</i> 8:7653. Liang SC, <i>et al.</i> 2003. <i>Eur. J. Immunol.</i> 33:27064. Barber DL, <i>et al.</i> 2006. <i>Nature</i> 439:6825. Keir ME, <i>et al.</i> 2005. <i>J. Immunol.</i> 175:73726. Koehn BH, <i>et al.</i> 2008. <i>J Immunol.</i> 181:5313
Gene ID	18566


Related Protocols

[Cell Surface Flow Cytometry Staining Protocol](#)


Other Formats

PE anti-mouse CD279 (PD-1), Purified anti-mouse CD279 (PD-1), PerCP/Cyanine5.5 anti-mouse CD279 (PD-1), APC anti-mouse CD279 (PD-1), Biotin anti-mouse CD279 (PD-1), FITC anti-mouse CD279 (PD-1), PE/Cyanine7 anti-mouse CD279 (PD-1), Brilliant Violet 421™ anti-mouse CD279 (PD-1), Brilliant Violet 605™ anti-mouse CD279 (PD-1), APC/Cyanine7 anti-mouse CD279 (PD-1), Brilliant Violet 785™ anti-mouse CD279 (PD-1), PE/Dazzle™ 594 anti-mouse CD279 (PD-1), Alexa Fluor® 647 anti-mouse CD279 (PD-1), Brilliant Violet 711™ anti-mouse CD279 (PD-1), GolnVivo™ Purified anti-mouse CD279 (PD-1), APC/Fire™ 750 anti-mouse CD279 (PD-1), Brilliant Violet 510™ anti-mouse CD279 (PD-1), Ultra-LEAF™ Purified anti-mouse CD279 (PD-1), APC/Fire™ 810 anti-mouse CD279 (PD-1) Antibody, PE/Fire™ 810 anti-mouse CD279 (PD-1) Antibody, PE/Cyanine5 anti-mouse CD279 (PD-1), PE/Fire™ 640 anti-mouse CD279 (PD-1)

Product Data



Con-A and IL-2 stimulated C57BL/6 mouse splenocytes (3 days) were stained with CD3 FITC and CD279 (clone 29F.1A12) APC/Fire™ 750 (top) or rat IgG2a, κ APC/Fire™ 750 isotype control (bottom).



For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587